

Coagulation Analyzer

Micro fluidic Coagulation analysis system (Electro chemical Method)

Fast:

Result delivered in 7 minutes

Easy Operation:

Single drop of whole blood, One-step rapid

Accurate Result:

Comparable diagnostic accuracy

Multiple parameters

✓ PT ✓ APTT

✓ FIB ✓ TT ✓ INR


*ACT Coming soon



HC100

Clinical application

Coagulation dysfunction screening



- Skin Petechiae, Purpura or Ecchymosis
- Gum bleeding
- Wound bleeding
- Gastrointestinal bleeding
- Mucosal bleeding
- Hematuria and Bloody stool
- Anemia

Pre-operative examination



- Patients with underlying Internal disease
- Pyometra
- Tumor
- Surgery
- Biopsy


Diagnosis and Prognosis of Internal Disease



- Multiple organ dysfunction and failure
- Obstetric diseases and complications
- Leukemia, sepsis and other blood diseases
- Immune system disease

Operating procedures

- 1 Click test
- 2 Scan QR code
- 3 Insert reagent strip
- 4 Add sample
- 5 Start detection



Pulse Diagnostics

Add: 4500 service road Burlington ontario Canada

Web: www.pulsediagnosics.net

Tel: 905333426

E-Mail: sales@pulsediagnosics.net